

### 3. Analysis of the effect of R-loops on replication fork progression

#### 0 day

- seed cells  $3 \times 10^5$  cells/well; 6-well plate; O/N
- to be at approx. 70% confluence next day;
- treatment with dox

#### 1st day

- pre-warm: DMEM, PBS
- cool PBS (in fridge or on ice)
- prepared all reagents and treatment mix in 15-mL falcons or 1.5mL-ependorf tubes, keep in water bath (37°C)

- CldU 3  $\mu$ L/mL + dox ( $\pm$ )
- incubate (30 min; 37°C; 6% CO<sub>2</sub>)
- wash pre-warm PBS (2 mL)

- pre-warm: Trp
- cool down centrifuge (4°C)

- IdU 12  $\mu$ L/mL + dox ( $\pm$ )
- + corresponding treatment
- incubate (30 min; 37°C; 6% CO<sub>2</sub>)
- wash pre-warm PBS (2 mL)

- bucket with ice; label 1.5mL-ependorf tubes


- harvest cells
  - 400  $\mu$ L Trp (approx.. 5 min; 37°C; 6% CO<sub>2</sub>)
  - add 1 mL DMEM to neutralize Trp
- pipet  $\uparrow\downarrow$  (3x times)
- transfer cells to 1.5mL-ependorf tube (labeled)

- spin (4 min; 4°C; 1500 rpm)
- aspirate supernatant
- resuspend 1 mL cold PBS, pipet  $\uparrow\downarrow$  (3x times)
- spin (4 min; 4°C; 1500 rpm)
- aspirate supernatant
- resuspend cells 500  $\mu$ L cold PBS  $\mu$ L;

- wipe glass slides with dH<sub>2</sub>O, label it with pencil
- CVs: 2 replicates / condition (well)
- prepare lysis buffer

- briefly vortex cells before counting
- count cells
- adjust volume to  $2.5 \times 10^5$  cells/mL

- 3  $\mu$ L of cells + 9  $\mu$ L lysis buffer
  - first the cell suspension – to approx. 2 CVs,
  - then add lysis buffer (\*mix ... start time)
- mix by gently with the tip (P200 pipette tip - yellow)
- avoid creation of bubbles
- lysis 5 min at r.t. ( starts after mixing);
- tilt the slide (15°- 45°)
- air dry the slides completely (15-20 min)
- place slides in to IF-cell(s)
- fix with fresh solution of methanol/acetic acid 3:1 (10 min; r.t.)
- dry O/N (keep the IF-cell open)



wash  
1. aspirate liquid  
2. add washing solution  
3. aspirate liquid  
..... repeat .....  
X. add desire liquid

treatment	c stock	c final	calculate	$\mu$ L / 1 mL
dox	1 $\mu$ g/mL	1 ng/mL		1 $\mu$ L
MUS81i	5 mM	5 $\mu$ M		1 $\mu$ L
CDDP	1 mM	10 $\mu$ M		10 $\mu$ L
CPT	100 $\mu$ M	100 nM		1 $\mu$ L
CldU	19 mM	57 $\mu$ M		3 $\mu$ L
IdU	28.24 mM	338 $\mu$ M		12 $\mu$ L

MM 1A (Master mix) in 15-mL falcon CldU 18  $\mu$ L / 6 mL

MM 2A in 15-mL falcon CldU 15  $\mu$ L / 5 mL

MM 1B in 15-mL falcon IdU 72  $\mu$ L / 6 mL

MM 2B in 15-mL falcon IdU 60  $\mu$ L / 5 mL

Cell culture room

Keep cells on ice




Trp: cells should look circular, be de-attached – floating; check progress under the microscope

? overly pipetting can cause lysis of cells

cold PBS (approx. 4°C)





CVs: First group: 10 CVs Second group: 8 CVs

IF-cell ... 75 mL



page 5

#### Counting cells



$$\text{cells/mL} = \frac{\# \text{ of cells counted}}{\# \text{ of big squares counted}} \times \text{dilution factor} \times 10\,000$$

$$V_1 c_1 = V_2 c_2$$

total volume = (calculated cells/mL x current volume) /  $2.5 \times 10^5$  cells/mL  
 volume to add = total volume – current volume in epp. 1/5

Lysis buffer				
solution	1 mL	solution	1 mL	2 mL
SDS 10%	50 $\mu$ L	SDS 20%	25 $\mu$ L	50 $\mu$ L
TrisHCl pH 7.4 -7.5	200 $\mu$ L	TrisHCl pH 7.4 -7.5	200 $\mu$ L	400 $\mu$ L
EDTA 500 mM	100 $\mu$ L	EDTA 500 mM	100 $\mu$ L	200 $\mu$ L
H <sub>2</sub> O sterile	650 $\mu$ L	H <sub>2</sub> O TK	675 $\mu$ L	1350 $\mu$ L

### 3. Analysis of the effect of R-loops on replication fork progression

#### 2<sup>nd</sup> day

- rehydrate slides in 1xPBS (2 x 5 min)
- denature DNA in 2.5 M HCl (1 h; r.t.)
- ✿ prepare and filter blocking solution
- ✿ calculate antibodies dilution
- + keep some filtered blocking solution for antibody dilution
- ✿ prepare 1L PBST 0.2%
- wash PBS (5 x 3 min)

#### Blocking

- blocking solution (recipe bellow)
- blocking in IF-cell (20 min)
- ✿ can be reuse multiple times (collect, keep at 4°C)

#### Primary antibody

- drain slides briefly dipping on a tissue
- dilute the antibodies in blocking solution (mix in the one ependorph)

anti-rat CldU 1:500  
 anti-mouse IdU 1:100  
 65 µL/CV; covered by cover slips  
 incubation (2 h; r.t.; the black staining box filled with H<sub>2</sub>O)

- wash in PBST 0.2% (5 x 3 min)
- blocking solution (leave while preparing the sc. Ab)

#### Secondary antibody

- dilute the antibody in blocking solution (mix in the one ependorph)

???  
 secondary antibody have to be kept in dark from on

anti-mouse Alexa-488 1:300  
 anti-rat cy3 1:150  
 65 µL/CV; covered by cover slips  
 incubation (1h; r.t.; the black staining box filled with H<sub>2</sub>O)

- wash in PBST 0.2% (5 x 3 min)
- 1xPBS (2 x)
- pour milliQ-H<sub>2</sub>O to the IF-cell and immediately take in out
- let slides air dry completely (15-20 min)
- mount with Fluoromont 30 µl/slide
- dry (approx. 1h)


#### Microscopy

- DM6000; oil 63x
- 7 fibers/image → 150 fibers .... approx. 40 images (20 + 20 from each replicates)

#### Analysis

- 150 fibers (150 green + 150 red = 300 measurements)


2.5 M HCl			
acids general: pour acid into water			
32% HCl → 25ml/100ml; 37% HCl → 20ml/100ml, 35% HCl → 22ml/100ml			
Blocking solution			
solution	70 mL	90 mL	140 mL
1xPBS	68.6 mL	88.2 mL	137.2 mL
BSA	1.4 g	1.8 g	2.8 g
Tween 20	70 µL	90 µL	140 µL
2% BSA, 0,1% Tween 20, 1xPBS, 0,22 µL filtered; can be reused			
PBST 0.2 %			
solution			
1xPBS	1 L		
Tween 20	2 mL		



#### Fixation

The preservation of biological tissues from decay. It terminates any ongoing biochemical reactions; may also increase mechanical strength or stability.


#### denaturation



#### DNA Denaturation

The compact structure of the double helix is maintained by the hydrogen bonds between base pairs and the van der Waals interactions between the stacked bases. Various agents (heat, strong alkalis, urea, and formamide) weaken such forces and promote the separation of the strands. The resulting unwound polynucleotide strands adopt a random arrangement.


#### blocking



#### Blocking

reduce non-specific background noise by covering exposed binding sites.

#### primary ant




#### A primary antibody

an immunoglobulin that binds directly to a specific target antigen (protein or biomolecule) to detect, quantify, or locate it within a sample. They are produced in various host species (e.g., mice, rabbits) and come as monoclonal (single epitope) or polyclonal (multiple epitopes) types.


anti-rat CldU & anti-mouse IdU  
 • both is antibody against BrdU  
 • each has different affinity to CldU or IdU

#### secondary ant



#### Secondary antibodies

essential research tools that bind to primary antibodies, allowing visualization or detection of target antigens. They are conjugated with labels (fluorescent, enzymatic, or biotin) for signal detection and enhance sensitivity by amplifying the signal, as multiple secondaries can bind one primary antibody.



CldU 555 nm  
 IdU 488 nm

replication – normal speed

replication – slowed down

Acquire images and analyse

### 3. Analysis of the effect of R-loops on replication fork progression

**3<sup>rd</sup> day**

ImageJ/Fiji

*Protocol in detail in Analysis\_Fiji\_ImageJ*

### 3. Analysis of the effect of R-loops on replication fork progression

#### Treatment

##### CPT

- camptothecin
- fast, the most common
- inhibition of Top 1 → it leaves nicks; ↑ gH2AX; run off; SSB → DSB during replication
- ↑ negative torsion stress → ↑ R-loops serve as obstacle for replication
- (CPT works similar as depletion of Top1)

##### CDDP

- cisplatin
- adducts + oxidative stress; resolved by translesion synthesis and repriming

##### MUS81i

- inhibitor of MUS81
- block transcription-replication conflict resolution
- MUS81 ... endonuclease cut DNA strand → release of supercoil → replication can continue

##### dox

- doxycyclin
- induction of RNase H1 wild type – overexpression
- for conflict resolution

#### Reagent & solutions

##### PBS

- Phosphate-Buffered Saline
- 10x PBS ... stock; 10x concentrated
- 1x PBS ... as working solution, can be written as PBS in protocol

##### PBST

- PBS + tween 20 (0.2%)
- tween 20 ... detergent
- to improve washes

##### DMEM

Dulbecco's Modified Eagle Medium

a widely used basal medium for supporting the growth of many different mammalian cells

##### Trp

- trypsin/EDTA (in this case)
- to detach adherent cells from culture surfaces (trypsinization) and dissociate tissues into single-cell suspensions. a protease, breaks down adhesion proteins
- EDTA, a chelating agent, removes calcium/magnesium to weaken cell-to-cell adhesion and boost enzyme activity + stabilization of Trp

##### BSA

- bovine serum albumin

##### Primary antibodies

- anti-rat CldU & anti-mouse IdU
- both is antibody against BrdU – with different affinity to CldU or IdU
- CldU, IdU, BrdU ... thymidine analogs





Figure 1. RNA/DNA hybrids and single rNMPs in DNA are substrates for RNase H2.

DNA is in blue, RNA is red line and letter R. The red oval represents RNA polymerase.





Source of images:

DOI:10.1101/gad.573310

DOI:10.1083/jcb.201104003https://annualreport.nichd.nih.gov/2013/crouch.html


### 3. Analysis of the effect of R-loops on replication fork progression

 Treatment schedule and scheme




seeding ± dox

First group



+ dox: 2.5 µL dox / 2.5 mL DMEM  
- dox: only DMEM


Second group



+ dox: 4 µL dox / 4 mL DMEM  
- dox: only DMEM


CldU ± dox

First group



CldU: 18 µL CldU / 6 mL DMEM  
+ dox: 2.5 µL dox / 2.5 mL CldU-DMEM  
- dox: 3.5 mL CldU-DMEM


Second group



CldU: 15 µL CldU / 5 mL DMEM  
+ dox: 3.5 µL dox / 3.5 mL CldU-DMEM  
- dox: 1.5 mL CldU-DMEM


IdU ± dox treatment

First group



IdU: 72 µL CldU / 6 mL DMEM  
+ dox: 2.5 µL dox / 2.5 mL IdU-DMEM  
- dox: 3.5 mL IdU-DMEM

Second group



CldU: 72 µL IdU / 6 mL DMEM  
+ dox: 3.5 µL dox / 3.5 mL IdU-DMEM  
- dox: 1.5 mL IdU-DMEM

treatment

treatment

